

Boiling Range Distribution of Crude Petroleum by Gas Chromatography

ASTM D5307

Configuration of an Agilent Technologies 7890A Series Gas Chromatograph with a flame ionization detector (FID).

Method covers the determination of water-free crude petroleum through 538°C (1000°F). Material boiling above 538°C is reported as residue. Method is applicable to whole crude samples that can be solubilized in a solvent to permit sampling by means of a microsyringe.

Repeatability Table RSD = Relative Standard Deviation SimDist Retention Time Mix Diluted 10:1 in Carbon Disulfide				
Decane (n-C10)			Eicosane (n-C20)	
Runs	Ret. Time	Area	Ret. Time	Area
1	5.437	2311.95068	13.288	2397.39893
2	5.437	2336.93433	13.287	2422.95386
3	5.436	2342.67187	13.287	2429.69629
4	5.436	2326.44702	13.286	2406.72827
5	5.436	2348.88135	13.286	2431.33569
6	5.435	2356.42334	13.286	2443.46045
7	5.436	2363.50464	13.285	2444.37085
Average	5.436	2340.97332	13.286	2425.13491
St. Dev.	0.00069	17.70794	0.00098	17.69156
% RSD	0.01269%	0.75643%	0.00735%	0.72951%
Triacontane (n-C30)			Tetracontane (n-C40)	
Runs	Ret. Time	Area	Ret. Time	Area
1	18.375	2389.72947	22.113	2356.52832
2	5.437	2410.94873	22.114	2380.63574
3	5.436	2416.81421	22.112	2386.19092
4	5.436	2395.22510	22.112	2369.17065
5	5.436	2419.11182	22.111	2393.84570
6	5.435	2431.78955	22.111	2404.91284
7	5.436	2433.19067	22.110	2407.65479
Average	18.374	2413.82994	22.112	2385.56271
St. Dev.	0.00076	16.68065	0.00135	18.56278
% RSD	0.00411%	0.69104%	0.00608%	0.77813%